

MagStrap[®] Instructions

1 Cleaning

Please clean your frames with the manufacturer recommended solution to clear any particulates and remove any oils before attaching the **MagStrap**.

2 Fitting

Rest the **MagStrap** flat on a hard surface and whilst securely holding your loupes, slide the silicone tubes over each arm of the frames.

3 Adjusting

Slide the toggle to the desired position to ensure a reliable and comfortable fit.





Follow instructions for use



Do not use with hearing aid



This device is not sterile



200mm safe distance from patients or operators with pacemakers



This device contains magnets which should not be used if you have a pacemaker, or with patients with pacemakers.

Model	MG2A
Known Contraindications	Silicone Allergy or Sensitivity
Material	Silicone Polymer and Neodymium (N52) Magnets

Please scan here to visit our support page and download your manual before use.




Speak to the specialists.

UK + 44 (0) 193 232 00 64 | info@bryantdental
AUS + 61 (2) 720 895 92 | conciergeau@bryantdental
www.bryantdental

UK office: The Barns, Hilltop Farm, Lyne Lane, Chertsey, UK, KT16 0AW
Australia office: Suite 18, 13U/175 Lower Gibbes St Roseville NSW 2069

MagStrap[®]

Manufacturer: Bryant Dental, The Barns, Hilltop Farm, Lyne Lane, Chertsey, KT16 0AW

 Dentistry, reimagined.