

# IGNIS

---

## Warranty

### LIMITED LIABILITY POLICY

Bryant Dental will repair any LED or wiring issue arising from a manufacturing defect for a period of 5 years, without charge.

Batteries and filters are covered for a period of 2 years.

The liability will not exceed the amount you originally paid for the product. The warranty does not cover defects arising from negligence, product modification, unauthorised repairs, accidental damage, theft, loss or acts of God.

A **BRYANT DENTAL** PRODUCT




## Warranty Safety Instructions

### CONDITIONS

- » Only use supplied wall charger, using another charger will damage the batteries and invalidate the warranty
- » Do not immerse any of the IGNIS light parts in water
- » If any problems occur with the unit return to Bryant Dental for repair, attempting to repair yourself will void the warranty
- » Do not change or modify any of the existing components without authorisation from Bryant Dental
- » Do not leave the light on inside a pocket – it will cause the unit to overheat
- » Do not leave the batteries on exposed skin
- » Do not charge unattended or leave in environments exceeding 40°C
- » Do not attach multiple batteries together
- » Batteries should be stored in a cool, dry place
- » If batteries behave abnormally immediately contact Bryant Dental

### Need help?

Contact our support team today at: **[www.bryant.dental](http://www.bryant.dental)**

   +44 (0)1932 320064 | +44 (0)7990 113723